

## Isolate Native Protein Complexes with an Improved Co-IP Kit, Validated for Lower Quantities of Antibody

### Capture Antigens and Interacting Proteins Using Only 10 – 50µg of Antibody

- Sufficient reagents to perform 50 reactions using 25µl of immobilized antibody support
- Shorter protocol
- Less antibody required
- Purify target proteins without contamination by the antibody in the eluate



**Special Offer** Quote **Special Offer Reference 2101** to receive 10% discount on Pierce Co-IP Kits, until 31 May 2010\*

### Ordering Information

Product Code	Description	List Price	Offer Price*
26149	<b>Pierce Co-Immunoprecipitation Kit</b> <i>Sufficient to perform 50 reactions using 25µl of immobilized antibody support.</i> Kit Contents: AminoLink Plus Coupling Resin, 2ml • Coupling Buffer (20X), 25ml • Sodium Cyanoborohydride Solution (5 M), 0.5ml • Quenching Buffer, 50ml • Wash Solution, 60ml • IP Lysis/Wash Buffer, 2 x 50ml • Modified Dulbecco's PBS (20X), 25ml • Conditioning Buffer (100X), 5ml • Elution Buffer, 50ml • Lane Marker Sample Buffer (5X), 5ml • Pierce Spin Columns, Screw Cap, 50 ea. • Microcentrifuge Collection Tubes, 2ml, 100 ea. • Microcentrifuge Collection Tubes, 1.5ml, 50 ea. • Pierce Control Agarose Resin, 2ml	€315.62	€284.06
23620	<b>Pierce c-Myc-Tag IP/Co-IP Kit</b> <i>Sufficient material to perform 25 IP/Co-IP assays.</i>	€408.82	€367.94
23625	<b>Pierce Mammalian c-Myc Tag IP/Co-IP Kit</b> <i>Sufficient material to perform 25 IP/Co-IP assays.</i>	€433.18	€389.86
21516	<b>Pierce GST Protein Interaction Pull-Down Kit</b> <i>Sufficient materials for conducting 25 pull-down assays using a GST-tagged fusion protein as the bait</i>	€327.27	€294.54
21277	<b>Pierce His Protein Interaction Pull-Down Kit</b> <i>Sufficient materials for conducting 25 pull-down assays using a His-tagged fusion protein as the bait</i>	€327.27	€294.54
23610	<b>Pierce HA-Tag IP/Co-IP Kit</b> <i>Sufficient material to perform 25 IP/Co-IP assays.</i>	€408.82	€367.94
23615	<b>Pierce Mammalian HA-Tag IP/Co-IP Kit</b> <i>Sufficient material to perform 25 IP/Co-IP assays.</i>	€433.18	€389.86
26147	<b>Pierce Crosslink IP Kit</b> <i>Sufficient reagents to perform at least 50 immunoprecipitation reactions</i>	€367.26	€330.53

\* Special Offer valid until 31 May 2010. Please quote the special offer reference and all relevant product codes when placing your order to avail of this offer. Failure to quote the special offer reference may result in list pricing being applied to your order. Discounts are not cumulative. All other special offers and discounts do not apply. A full list of kit contents can be found at [www.medical-supply.ie/category\\_6.shtml](http://www.medical-supply.ie/category_6.shtml)