


TransFast™ Transfection Reagent


The Quick and Easy Solution

Features

- **Fast:** Transfect in 1 hour. Transfection times can be decreased to as little as 30 minutes with certain cell lines.
- **Easy to Use:** Resuspend the reagent in water, freeze, thaw, mix with DNA, and add to cells.
- **Efficient:** High-efficiency transfection – transient and stable – in many cells.
- **Robust:** Requires less optimization than other systems. Allows transfection of cell types such as primary cell cultures that require continuous exposure to serum.



Cell density and transfection efficacy of TransFast™ Transfection Reagent.



Expression of Monster Green™ Fluorescent Protein in transfected human neuroblastoma SH-SY5Y cells

To apply for your free sample, visit
www.medical-supply.ie/request_sample.shtml

Special Offer Quote **Special Offer Reference 2104** to receive discounted pricing until 31 May 2010*

Ordering Information

Product Code	Description	No of Transfections	List Price	Offer Price*
E2431	TransFast™ Transfection Reagent, 1.2mg	400	€345.00	€245.00

* Special Offer valid until 31 May 2010. Please quote the special offer reference and all relevant product codes when placing your order to avail of the special pricing. Failure to quote the special offer reference may result in list pricing being applied to your order. Discounts are not cumulative. All other special offers and discounts do not apply.