

# ZOOM

sebia

“Testing Simplified”

## Hydrigel 5 von Willebrand multimers

## Von Willebrand Factor multimers assay simplified




### INTENDED USE

Detection of the abnormalities in the VWF multimers distribution related to a hemorrhagic diathesis, as seen in types 2A and 2B VWD (von Willebrand disease) and AVWS (Acquired von Willebrand syndrome).



1 day



### MAIN BENEFITS

- Same day results
- Simplified method
- Clear and reproducible profiles
- Reduced hands-on time
- Ready to use gels and reagents
- All multimer sizes

### References

- PN 4359 Hydrigel 5 von Willebrand multimers kit
- 11 wells gel\*\*
- PN 4747 5 von Willebrand multimers visualization kit
- PN 1651 Hydrasys von Willebrand multimers accessories
- PN 1652 Metal Weight

### Availability


- All countries
- For Research Use Only\*\*\*

### Posters\*\*\*\*

- IFCC 2015
- EHA 2015
- ISLH 2016
- WFH 2016

### Instrument

- Hydrasys 2 + Gel scan
- Hydrasys 2 scan



### DENSITOMETRY

- Makes interpretation easy
- Allows the overlay with a normal control
- Enables estimating the relative concentration of each multimer subset

\* Normal control (not provided by Sebia)

\*\* Under development

\*\*\* Contact your Sebia representative for more details

\*\*\*\* Available on our website